



Braun Swing-A-Way®







Solves the problem of tight parking spaces with parallel boarding and its unique swivelling platform.







For your safety, the hydraulic roll stop raises to the vertical position before the platform leaves the ground.

The Braun Swing-A-Way creates a unique ability to park close to buildings or narrow pavements as the wheelchair user rolls onto the platform in parallel to the vehicle. The telescopic platform extends to the necessary length as the wheelchair user rolls onto it.

The lifts are fully automatic in operation and operated by a user or care-giver using either the standard handheld control or the controls located on the lift.

Quick facts

- Requires a minimal parking space. The deployed platform protrudes no more from the side of the vehicle than the fully opened front door.
- Does not obstruct normal use of door opening
- Platform handrail provides added security for wheelchair user
- EMC test approved
- Type approved
- CE-marked

To guarantee your safety all Autoadapt products are designed and tested according to current directives and standards.

Specifications

CONSTRUCTION Steel with powder coat finish

OPERATING TEMPERATURE -23°C to +65°C









Braun Swing-A-Way®

Dimensions


Side view




Side view



Top view



Top view



		Pomer suf	Clear platfor	in Clear platfor	r Overallik	Mat ground	trat ground overallift	thin openio	oor din dear d		n Lifting capacity	ited Unit weight
Swing-A-Way			Α	В	С	D	Е	F	G	Н		
700115	Swing-A-Way L800XEVA	12 V	66	117	117	69	115	104	118	71	363	122
700117	Swing-A-Way L800XVA	12 V	74	92	117	69	117	104	118	71	363	111